ParaPost®XSystem

ParaPost®XP

ParaPost®Plus

ParaPost®

One-Office-Visit Technique


ParaPost®XH

One-Office-Visit Technique

ParaPost®XP

ParaPost®

In-Direct Casting Technique

Instructions for Use

CE 0344

 **COLTENE**

EN	Instructions For Use	3
DE	Gebrauchsanweisung	12
FR	Mode d'emploi	20
ES	Instrucciones de uso.....	27
IT	Istruzioni per l'uso	35
NL	Gebruiksaanwijzing	42
SV	Användarinstruktioner.....	50
DA	Brugsanvisning	57
FI	Käyttöohjeet.....	64
NO	Bruksanvisninger	71
PT	Instruções para uso	78
EL	Οδηγίες χρήσης	86

ParaPostXP, ParaPost Plus, ParaPost

One-Office-Visit Technique

Intended Use: Metal posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

 CAUTION!	Federal law (United States) restricts this device to sale by or on the order of a healthcare practitioner.
 WARNING!	Failure to properly follow the directions in this IFU may result in patient injury or compromise the integrity of the restoration.
 NON STERILE	Posts are delivered non-sterile and are to be sterilized prior to use.
 WARNING!	Device should not be used in patients with known material sensitivity.

Note: Posts are single use devices and should never be salvaged for re-use. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture. For illustrations see pages 100 - 102.

Note that titanium posts have virtually the same radiopacity as composite and gutta-percha. Careful scrutiny of the radiograph will allow for the differentiation of post from gutta percha and composite.

We urge the use of Hygenic® Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 01**).^{*} If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.^{*}
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)^{*} to prevent perforation or weakening of root walls, and at least 4-5 mm of gutta-percha must be left intact to protect the apical seal.
3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat with Hygenic Gutta-Percha before preparing the post space.

4. To begin paralleling the post space, select the ParaPost Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 01**). (See "Approximate Diameter Comparisons" chart, **Fig. A**). To help maintain the preplanned depth of the preparation, ParaPost Drills have been pre-marked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contraangle (**750 - 1,000 RPM**). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Drill until the preplanned diameter and depth are achieved (**Fig. 02**).
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (**Fig. 03**).
7. Select the ParaPost post that corresponds to the last drill used to prepare the post space. Remove head portion from the post. The head portion is for identification purposes only. Insert post into post space (**Fig 04**). Check for occlusal clearance. Remove the post; shorten if necessary.
8. To prepare for cementation, clean and rinse the canal (**Fig. 05**) and dry with Hygenic Paper Points. Use the dental cement of your

choice according to the manufacturer's instructions.

9. Use a Lentulo spiral to coat the walls of the post space with cement (**Fig. 06**).
10. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent (**Fig.07**).
11. If necessary, to obtain additional occlusal clearance after the cement has set, make any adjustments needed to the post head while holding the post with a hemostat (**Fig. 08**).
12. After the cement has set, remove any excess before completing the core and fabricating the final restoration (**Fig. 09**).

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Alternatively, sterilize in a pre-vacuum chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130 degrees C for 15 minutes, following standard procedures. Alternatively, sterilize in a pre-vacuum chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.


Warning Notice: DO NOT BEND POSTS.

**References and/or literature available upon request.*

ParaPostXH

One-Office-Visit Technique

Intended Use: Metal posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

Note: Posts are single use devices and should never be salvaged for

re-use. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture. For illustrations see page 100 - 102

Note that titanium posts have virtually the same radiopacity as composite and gutta-percha. Careful scrutiny of the radiograph will allow for the differentiation of post from gutta percha and composite.

We urge the use of Hygienic® Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 10**)* If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.*
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)* to prevent perforation or weakening of root walls, and at least 4-5 mm of guttapercha must be left intact to protect the apical seal
3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat with Hygienic Gutta-Percha before preparing the post space.

- To begin paralleling the post space, select the ParaPostX Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (Fig. 10). (See "Approximate Diameter Comparisons" chart, Fig. A). To help maintain the preplanned depth of the preparation, ParaPostX Drills have been pre-marked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPostX Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle (750 - 1,000 RPM). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

- Sequentially step up to the next larger ParaPostX Drill until the preplanned diameter and depth are achieved (Fig. 11).
- Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (Fig. 12).
- Select the ParaPostXH post that corresponds to the last drill used to prepare the post space. Insert post into post space (Fig. 13). Check for occlusal clearance. Remove post; if post has to be shortened, make adjustments from apical end (Fig. 14). Reamfer the end of the post to its original shape.
- To prepare for cementation, clean and rinse the canal (Fig. 15) and dry with Hygienic paper points. Use the dental cement of your

choice according to the manufacturer's instructions.

- Use a Lentulo spiral to coat the walls of the post space with cement (Fig. 16).
- Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent (Fig. 17).
- If necessary, to obtain additional occlusal clearance after the cement has set, make any adjustments needed to the post head while holding the post with a hemostat (Fig. 18).
- After the cement has set, remove any excess before completing the core and fabricating the final restoration (Fig. 19).

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Alternatively, sterilize in a pre-vacuum chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130 degrees C for 15 minutes, following standard procedures. Alternatively, sterilize in a pre-vacuum chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning Notice: DO NOT BEND POSTS.

**References and/or literature available upon request.*

ParaPostXP, ParaPost

Casting Technique - Indirect

Note: Posts are single use devices and should never be salvaged for re-use. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture.

For illustrations see page 100 - 102

We urge the use of Hygienic® Rubber dam with these procedures.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 20**).^{*} If that amount of tooth structure is not available,

it should be created by periodontal crown lengthening and/or orthodontic extrusion.^{*}

2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)^{*} to prevent perforation or weakening of root walls, and at least 4-5 mm of guttapercha must be left intact to protect the apical seal.
3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat with Hygienic Gutta-Percha before preparing the post space.
4. To begin paralleling the post space, select the ParaPost Drill which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 20**). (See "Approximate Diameter Comparisons" chart, Fig. A). To help maintain the preplanned depth of the post space, ParaPost Drills have been premarked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contraangle (**750 - 1,000 RPM**). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the possibility of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Drill until the previously determined diameter and depth are achieved (**Fig. 21**).
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box without undercuts (**Fig. 22**).

Indirect Only

7. Place into the prepared post space the ParaPost Impression Post corresponding to the last drill used. If necessary, adjust the length of the Impression Post so the post does not contact the impression tray.
8. Take an impression with an elastomeric impression material (i.e., President). Make sure the anti-rotational box is filled with impression material

NOTE: To ensure retention of the post in the impression material when using Hydrocolloid, form a nail-head on the top of the Impression Post by pressing on the head with a hot instrument. DO NOT ALTER the plastic Impression Post when using silicone or rubber-based impression material.

9. Insert the corresponding ParaPost Temporary Post into the prepared post space. If necessary, shorten from the apical end. Fabricate a provisional crown. To preserve the integrity of the post space, apply temporary cement ONLY to the margins of the crown; DO NOT PLACE any cement into the post space or on the post.
10. Pour a model and fabricate a precision removable die using the Pindex® or Accutrac® System.

11. Place the corresponding ParaPost Burnout Post Into the post space of the model. If necessary, shorten the Burnout Post from the coronal end.

12. Complete the core pattern.

13. Carefully remove the post/core pattern

14. Sprue, invest and cast. (When using a plastic ParaPost Burnout Post, cast with any alloy)

15. To prepare for cementation, clean and rinse the canal and dry with Hygenic Paper Points. Use the dental cement of your choice according to the manufacturer's instructions.

16. Use a Lentulo spiral to coat the walls of the post space with cement.

17. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent.

18. In order to obtain maximum retention to crown, DO NOT POLISH the casting. After the cement has set, remove any excess before completing the core and fabricate the final restoration.

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Alternatively, sterilize in a pre-vacuum

chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Temporary Posts: Prior to each use, sterilize by autoclave, following standard procedures.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning Notice: DO NOT BEND POSTS.

**Reference and/or literature available upon request.*

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Coltène/Whaledent Inc.
Coltène/Whaledent AG
Coltène/Whaledent GmbH + Co. KG