Carefully read instructions prior to use

**Product description:**

**GrandioSO Flow** is a flowable composite. In contrast to conventional flows, it is characterized by its increased filler content and stability. The radiopaque **GrandioSO Flow** contains 81% w/w (**GrandioSO Flow WO** contains 82% w/w) inorganic fillers in a methacrylate matrix and cures under blue light. In **GrandioSO Flow WO** an extremely radiopaque version is available. **GrandioSO Flow** is used with a dentine-enamel bond.

**Shades:**

A1, A2, A3, A3.5, A4, B1, B2, C2, I, OA2, BL, WO, GA5

**Indications:**

- filling minimally invasive cavities of all classes
- filling small class I cavities and extended fissure sealing
- filling class II – V cavities including V-shaped defects and cervical caries
- blocking out undercuts
- lining or coating cavities
- repairing fillings and veneers
- luting translucent prosthetic pieces (e.g., full ceramic crowns, etc.)

**Contraindications:**

**GrandioSO Flow** should not be used if the prescribed application technique cannot be employed or if there is a known allergy to any of its ingredients (methacrylates, BHT).

**Preparation:**

Clean the teeth to be treated.

**Shade selection:**

To select a shade, compare the enclosed shade guide with the cleaned tooth before anesthesia. Moisten the shade sample and tooth and select the shade, preferably in daylight.

**Creating a dry working field:**

Ensure that the working field is sufficiently dry. Use of a rubber dam is recommended.

**Cavity preparation:**

In principle, cavity preparation should be carried out according to the rules of adhesive filling therapy and should be minimally invasive to conserve healthy tooth substance. Bevel the enamel margins on anteriors. Pay attention to special preparation shapes for extended fissure sealing and deciduous teeth!  
Clean and dry the cavity afterwards. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

**Lining:**

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

**Matrix placement:**

Take appropriate measures for areas where the applied **GrandioSO Flow** can flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

**Handling:**

Allow the material to reach room temperature before use.

**Bonding material:**

**GrandioSO Flow** is used in the adhesive technique with a dentine-enamel bond. All light-curing bond materials can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

**Application of GrandioSO Flow:**

**Syringe:** Lock the enclosed application tip on the syringe with a clockwise turn and directly apply **GrandioSO Flow**. For fillings with a thickness of over 2mm, apply and cure in layers. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe.

**Caps:** Remove protective cap. Insert Caps into the opening of the dispenser following the respective instructions for use. Turn Caps to the desired direction and remove the protective cover. Apply **GrandioSO Flow** directly into the cavity by slowly and evenly pressing the levers of the handle together.

**Light-curing:**

Conventional polymerization devices are suitable for light-curing the material. The curing time is as follows when using an LED/halogen light with a minimum light-output of 500 mW/cm<sup>2</sup>.

20 s for shades A1, A2, A3, A3.5, B1, B2, C2, BL, I  
40 s for shades A4, OA2, WO, GA5

Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced. Incomplete curing may lead to discoloration and discomfort.

**Finishing:**

The restoration can be finished and polished immediately after removing the matrices (e.g. fine or extra-fine diamond grinder, polisher). The tooth should be fluoridated as a final step.

**Information/precautionary measures:**

- No known side effects. Sensitization in hypersensitive patients, however, cannot be excluded.
- The enclosed application tips are intended for a single use and must be disposed of after only one use. Return the cap to the syringe after use.
- Phenolic substances, especially preparations containing eugenol or thymol, interfere with the curing of filling composites. Avoid the use of zinc oxide eugenol cements or other eugenol-containing materials in combination with filling composites.

**Storage:**

Store at temperatures between 39 °F - 73 °F (4 °C - 23 °C). To avoid exposure to light and potentially resulting polymerization, syringes should be closed immediately after dispensing. Do not use **GrandioSO Flow** after expiry date.

**PRODUCT ORDERING INFORMATION:**

**GrandioSO Flow**

Kit 5 x 2.0gm syringes (A1, A2, A3, A3.5, WO)

REF 2720

Pack 2 x 2.0gm syringes in shades

A1	REF 2722	C2	REF 2732
A2	REF 2723	I	REF 2733
A3	REF 2724	OA2	REF 2734
A3.5	REF 2725	BL	REF 2736
A4	REF 2726	WO	REF 2737
B1	REF 2729	GA5	REF 2727
B2	REF 2730		

**GrandioSO Flow Caps**

Pack 16 x 0.25gm caps in shades

A1	REF 2742	A3.5	REF 2745
A2	REF 2743	A4	REF 2746
A3	REF 2744	B1	REF 2749

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.  
For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VOCO GmbH Phone +49 (4721) 719-70