

THE DENTAL
SOLUTIONS
COMPANY™



2021 PRODUCT CATALOG

Pioneers in Performance



[dentsplysirona.com](https://www.dentsplysirona.com)

Product Returns Terms & Conditions:

All returns must receive a RETURN MATERIALS AUTHORIZATION (RMA) number from the Dentsply Sirona Customer Service Department prior to shipment. This RMA number must be clearly noted on the outside of the shipping carton. If products are shipped without an RMA, the customer will be charged a processing fee.

Terms and Conditions apply. Product availability varies by market. Visit us online at dentsplysirona.com or contact your local customer care for more details.

Visit us online at dentsplysirona.com where you can shop our full line of products, download how-to-guides, get information on patient care and clinical education, watch videos and more.

Products



SureSmile® Aligner

03-05

SureSmile Digital Workflow	4
Why SureSmile?	5



Plastic

06-12


Plastic Selection Chart	7
Essix PLUS: Retainers and Bruxism Appliances	8
Essix ACE: Aligners and Retainers	9
Essix C+: Aligners and Retainers for Bruxers	10
Essix Nightguard Laminate: Bruxism Appliances.....	10
Essix TRAY RITE: Whitening Trays	10
Essix A+: Bite Planes.....	11
Essix Dual Laminate: TMJ Splints	11
Essix Pontics: Temporary Bridges.....	12
Essix Laminated Sports Mouthguard Material: Contact Sports.....	12



Thermoforming Equipment

13-14

Drufomat Scan Pressure Machine.....	14
Pressure Polymerization Insert, Granules, Funky Tool.....	14



Trimming Tools

15-20

Laboratory Discs	16
Precision Electric Heat Knife.....	17
Sports Mouthguard Trimming Bur	17
Finishing Liquid.....	17
Appliance Trimming Kit.....	18
Trim-Rite Release Agent.....	18
Accentuator Tool.....	18
Freeze Spray Coolant.....	18
Cordless Dremel	19
Essix Lab Knife	19
Scissors and Shears	19
Polishing Brushes & Border Polishers.....	20
Burs & Trimming Wheel.....	20
Micrometer.....	20



Interproximal Reduction

21-26

SpaceFile Contra-Angle Handpiece & Air Motor	22
SpaceFile IPR.....	23
Flexview, Diamond, & Serrated IPR Strips	24
Qwikstrip IPR	25
Autoclavable Aluminum Tube	25
IPR Gauge Set.....	25
IPR Burs.....	26




Patient Aids

27-29

Retainer Brite Cleaning Tablets	28
Sonic Cleaner.....	28
Chewies Aligner Tray Seaters.....	29
Appliance Remover Tools	29
Therapy Wafers	29



SureSmile® Aligner



THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

SureSmile® Digital Workflow

- Simple, easy-to-use digital process
- Results in biologically achievable targets and 50% less refinements than competitive aligners*
- Web-based software, accessible from any device
- Accepts STL files from most intraoral scanners

A Complete System Designed for Ultimate Versatility

The SureSmile Digital Workflow puts you in control.

Choose your treatment plan:

Create your own consult with SureSmile Digital Lab experts

Choose your fabrication:


SureSmile Digital Lab or make in-house

Choose your pricing plan:

Flat-rate, unlimited plan or pay-as-you-go



SCAN
Primescan™



PLAN
SureSmile Software



FABRICATE
SureSmile Arch Models



TREAT
SureSmile Aligners

*based on data analysis of SureSmile® Aligner global refinement rate compared to data reported in the survey of Klein/Vogels/Vogels, 2020 JCO Study of Orthodontic Diagnosis and Treatment Procedures (JCO/October 2020, p. 581 et seq.).



Need a digital scanner?

Discover **Primescan™**


As part of the Dentsply Sirona family, Primescan is a preferred scan source that offers best-in-class digital technology. The speed, accuracy and ease-of-use Primescan offers is unmatched. Plus, as is consistent with SureSmile Aligner's commitment to training and support, Primescan is backed by world-class educational resources and a user community to help guide you every step of the digital way.

Experience the simplicity of scanning.

Contact your Dentsply Sirona rep or learn more at dentsplysirona.com/primescan

Why SureSmile?

Combining clinically advanced planning capabilities with state-of-the-art manufacturing, SureSmile® planning software delivers a complete, clinician controlled, clear aligner treatment solution for you and your patients.



Freedom of choice

User-friendly, open architecture platform accepts files from most intraoral scanners



Versatile software

Robust software offering clinicians options to choose their desired level of treatment control; from single aligner to comprehensive, complex orthodontic cases.



Clinical support

Clinicians have confidence in SureSmile® knowing they have a team of dedicated specialists who support them in the Digital Lab and online through SureSmileU.


Support from the world's most advanced digital orthodontics platform.

Our support services combine the clinical expertise and technical skill of SureSmile orthodontists and technicians to provide clinically advanced digital modeling and appliance planning as you need it. You'll receive support from our expert digital lab service team as well as from the extensive training resources online. Plus, we're just a phone call away whenever you need assistance.





Plastic



Appearance. Strength. Durability.

There are plenty of reasons to insist on authentic Essix® Plastics. But is there a better reason than because the material you use is an integral part of your treatment plan?

Dentsply Sirona is a highly-respected name in dental plastics and polymer technology. When it comes to meeting the needs of clinician and patient alike, Essix® Plastics give you the quality you demand and the versatility you need.

Choosing the Right Plastic

	Aligners – Standard	Aligners – Bruxers	Retainers – Standard	Retainers – Bruxers	Whitening Trays	Nightguards	Mouthguards – Contact Sports	Bruxism Appliances	Temporary Bridges	TMJ Splints	Accepts Acrylics/ Bonding Agents	Bite Planes
Essix® PLUS™ Plastic			X	X				X	X			
Essix® ACE® Plastic	X		X						X			
Essix® C+® Plastic		X		X								
Essix® Nightguard Laminate						X						
Essix® Dual Laminate										X	X	
Essix® A+® Plastic											X	X
Essix® TRAY RITE® Plastic					X							
Essix® Laminated Sports Mouthguard Material							X					

Chart indicates primary use.

Essix® plastics are not made with natural rubber latex or BPA.

The Original Essix® Clear Retainer


Never outperformed.

ESSIX® PLUS™ Plastic is proven to look clearer, last longer and resist more staining than leading competitors. That means greater patient compliance and better case outcomes. Its proprietary formulation is user friendly. Appliances remove easily from models and trim effortlessly.

43% more clear than leading competitors*

Clarity

on a thermoformed sheet




*Data on file.

Up to 2X more durable than leading competitors*

Durability


Energy to Break (Joules)



Up to 7X more stain-resistant than leading competitors*

Staining



n Before Soak n After Soak



Up to 51% more wear-resistant than leading competitors*

Wear-Resistance

Material Loss in mg



Primary Use: Retainers and Bruxism Appliances

ESSIX® PLUS™ Plastic

- Excellent clarity and durability
- Easy to form and trim, excellent wear and impact resistance
- Accepts acrylics/bonding agents
- Protective film on both sides
- Also use for aligners and temporary bridges
- Patent #5,692,894


THICKNESS	SHAPE	ITEM #	QTY
.035" 7/8mm	5" Square	1APL35S5	Box of 80
	125mm Square	1APL35S125	
	125mm Circle	1APL35C125	
	120mm Circle	1APL35C120	
.040" 1mm	5" Square	1APL40S5	Box of 80
	125mm Square	1APL40S125	
	125mm Circle	1APL40C125	
	120mm Circle	1APL40C120	

Essix® plastics are not made with natural rubber latex or BPA.

Essix® ACE® Performance Plastic for Aligner Treatment

Proof is in the Plastic.

Essix ACE has been straightening teeth for over 20 years, and is proven to be the clear choice for aligner plastic.



Primary Use: Aligners and Retainers


ESSIX® ACE® Plastic

- Clear material has protective film on both sides
- Store in a cool, dry place
- Also use for temporary anterior bridges
- Patent #5,692,894

Essix ACE .030¹

Maintains its gentle forces
UP TO **14 DAYS**

Essix ACE can hold its forces for longer than 14 days during treatment, if needed.



THICKNESS	SHAPE	ITEM #	QTY
.030" 3/4mm	5" Square	1A350	Box of 100
	125mm Square	1A347	
	125mm Circle	1A325	
	120mm Circle	1A320	
.035" 7/8mm	5" Square	1AH50	Box of 100
	125mm Square	1AH347	
	125mm Circle	1AH25	
	120mm Circle	1AH20	
.040" 1mm	5" Square	1A450	Box of 100
	125mm Square	1A478	
	125mm Circle	1A425	
	120mm Circle	1A420	

1. Essix ACE .030¹ consistently applies more force than a multi-laminate plastic.

Essix® ACE® Advantages



COST EFFECTIVE

Less expensive per sheet than a multi-laminate plastic.



TREATS A WIDE RANGE OF MALOCCLUSIONS

...due to more available force.



HAPPIER PATIENTS

Can result in less office visits, shorter treatment times and improved patient comfort.*



BIO-COMPATIBLE


Safe for the patient and clinician for its intended use.

*For In-office aligner design adding more aligners to the treatment plan or creating smaller increments of tooth movements in the treatment plan will improve patient comfort.

Primary Use: Aligners and Retainers for Bruxers

ESSIX® C+® Plastic

- Contact clarity
- Store in a cool, dry place
- Not for use with acrylics/bonding agents
- Also use for bite planes
- Patent #5,692,894




THICKNESS	SHAPE	ITEM #	QTY
.040" 1mm	5" Square	1C500	Box of 100
	125mm Square	1C478	
	125mm Circle	1C125	
	120mm Circle	1C120	

Primary Use: Bruxism Appliances

ESSIX® Nightguard Laminate

- Single-sided construction; soft surface
- Excellent shape memory and abrasion resistance
- Also use for orthodontic repositioners and periodontal splints



THICKNESS	SHAPE	ITEM #	QTY
.120" 3mm	5" Square	1MENL	Box of 12

Primary Use: Whitening Trays

For Use With Vacuums


ESSIX® TRAY RITE® Plastic

- Clings to the tooth surface
- Doesn't stretch like EVA material - helps prevent bleaching material from leaking into mouth
- No reservoirs needed when fabricated on a vacuum machine
- Patent #6,089,869



THICKNESS	SHAPE	ITEM #	QTY
.020" 1/2mm	5" Square	70500	Box of 50
	125mm Circle	70125	
	120mm Circle	70120	

Primary Use: Bite Planes



ESSIX® A+® Plastic

- Clear material has protective film on both sides
- Accepts acrylics/bonding agents
- Store in a cool, dry place
- Additional uses include splints, nightguards (2mm and 3mm) and snoring appliances
- Patent #5,692,894

THICKNESS	SHAPE	ITEM #	QTY
.020" 1/2mm	5" Square	10025	Box of 100
	125mm Square	10020	
	125mm Circle	11020	
	120mm Circle	11021	
.030" 3/4mm	5" Square	10035	Box of 100
	125mm Square	10030	
	125mm Circle	11030	
	120mm Circle	11031	
.040" 1mm	5" Square	10045	Box of 100
	125mm Square	10040	
	125mm Circle	11040	
	120mm Circle	11041	
.060" 1 1/2mm	5" Square	10065	Box of 50
	125mm Square	10060	
	125mm Circle	11060	
	120mm Circle	11061	
.080" 2mm	5" Square	10085	Box of 50
	125mm Square	10080	
	125mm Circle	11080	
	120mm Circle	11081	
.120" 3mm	5" Square	13500	Box of 25
	125mm Circle	13125	
	120mm Circle	13120	


Primary Use: TMJ Splints



ESSIX® Dual Laminate


- Dual-sided construction: hard top surface withstands mastication pressure while soft bottom surface provides patient comfort
- Accepts bonding agents
- Also use for nightguards, clenchers and grinders
- Note: Not recommended for use in a Drufomat

THICKNESS	SHAPE	ITEM #	QTY
.080" 2mm	5" Square	12500	Box of 12
	125mm Circle	12125	
.120" 3mm	5" Square	1MEDL	Box of 12
	125mm Circle	1MEDL125	



ESSIX® Pontics

- Wear-resistant, vacuum-processed synthetic resin
- Special lingual channel ensures pontic locks into the plastic appliance
- Patent #8,070,485



* Vita A, B & C designations correspond to the Vita® shade guide, which is meant to be a guide, not a match. Vita is not a trademark of Dentsply Sirona.

Upper Laterals EUT6-(color)

Mold Classification	Tapered	
Length Lateral (B)	9.9mm	
Width Lateral (C)	6.2mm	

Tapered Upper Anteriors TUA6-(color)

Mold Classification	Tapered	
Total Width on Curve (A)	49.5mm	
Length Lateral (B)	10.8mm	
Width Lateral (C)	8.3mm	

Square Upper Anteriors SUA6-(color)

Mold Classification	Square	
Total Width on Curve (A)	57mm	
Length Central (B)	11mm	
Width Central (C)	9.4mm	

Lower Anteriors ELA6-(color)

Total Width on Curve (A)	37.5mm	
Length Central (B)	9.7mm	
Width Central (C)	5.1mm	

0° Posteriors - Upper	
PUMO-(color)	
Width Upper (D)	32.2mm

0° Posteriors - Lower	
PLMO-(color)	
Width Lower (E)	34mm

ESSIX® PONTICS SET **EPO04-W or EPO04-Y**
Includes 6 cards, 1 ea. type



4mm only


4mm & 5mm

ESSIX® Laminated Sports Mouthguard Material

- Two laminated layers of protection
- A 1/8" plastic insert helps protect the incisal edges of teeth
- Store in a cool, dry place
- Pre-laminated

SHAPE	THICKNESS	COLOR	ITEM #	QTY
5" Square	.160" 4mm	CLEAR:	1MLS5	Box of 12
		1-COLOR:	1MLS5-(color)	
		2-COLOR:	1MLD5-(color)	
	3-COLOR:	1MLT5 Red/White/Blue		
	.200" 5mm	CLEAR:	1MLS55	Box of 12
1-COLOR:		1MLS55-(color)		

Essix® plastics are not made with natural rubber latex or BPA.



Thermoforming Equipment



THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

Drufomat Scan Pressure Machine

- Fabricate top-quality nightguards, sports mouthguards, retainers, and more
- Easy operation: barcode scanner programs heating and cool-down times
- Variable pressure of up to 6 bars (87 psi) produces excellent internal adaptation
- Vertical forming design ensures consistent results
- Dimensions: 18.1" h x 11.4" w x 11.4" d
- Two-year warranty

3300 Drufomat Scan, 115v, Red
3300S Drufomat Scan, 115v, Silver

Free With Machine

- Guidebook
- Quick Connection Kit
- Plastics Guide
- Plastic Sample Kit Includes:
 Essix ACE® Plastic, Essix C+® Plastic,
 Essix® Sports Mouthguard Material



Drufomat is not a trademark of Dentsply Sirona.



Optional Pressure Polymerization Insert

- Pressure pot for cast flasks (without water)
- For repairs or extensions of dentures (with water)


3227



Optional Granules

- Allows use of models with higher or untrimmed bases when thermoforming
- Two bottles of granules are required to fill the F insert reservoir

4282 1.2kg Bottle, Lead



Funky Tool

- Use the Funky Tool with the Drufomat Scan to create custom plastic sheets
- Mix and match colors to create totally awesome patterns
- Let your patients choose their favorite team colors
- There are all sorts of funky things you can do!
- Six-piece set includes Cutting Table, Cutting Guide, Heating Pan, Chimney, Unit Adapter, and Pressing Lid

3229 Funky Tool




Trimming Tools



Lab Essentials

All discs are for straight handpieces only and include a built-in mandrel unless otherwise noted.

Autoclavable; do not immerse in cold-sterilization solution.



Diamond Discs

- Laboratory disc for separating and contouring of ceramics
- Built-in mandrel
- For straight handpiece only

	911HHP140	911HHP220
Size	14mm	22mm
Sided	Double	Double
Coverage (B)	2mm	3mm
Thickness (L)	0.15mm	0.15mm
Max Speed	30,000 RPM	20,000 RPM

Flexview™ Mesh Discs

- Laboratory disc for ultra-fine contouring/shaping on plaster, ceramics and acrylics
- Not for breaking contact
- Built-in mandrel
- For straight handpiece only

	934HP140	934HP180	934HP220
Size	14mm	18mm	22mm
Coverage (B)	2mm	3mm	3mm
Thickness (L)	0.18mm	0.18mm	0.18mm
Max Speed	30,000 RPM	25,000 RPM	20,000 RPM

Flexview™ Curved Discs

- Laboratory disc for superfine separating and contouring of ceramics
- Curved slots provide visibility
- Built-in mandrel
- For straight handpiece only

	983HP220
Size	22mm
Sided	Double
Coverage (B)	7.0mm
Thickness (L)	0.10mm
Max Speed	20,000 RPM

Flexible Serrated Discs

- Laboratory disc for separating and contouring acrylics
- Built-in mandrel
- For straight handpiece only


	946HP180	946HP220
Size	18mm	22mm
Sided	Double	Double
Coverage (B)	3mm	3mm
Thickness (L)	0.20mm	0.20mm
Max Speed	25,000 RPM	20,000 RPM

Tools to create custom sports mouthguards.

Precision Electric Heat Knife

- Best for trimming custom sports mouthguards and scalloping custom bleach trays
- Melt and cut thermoforming plastic – up to 4mm splint material
- Built-in regulator
- Three-month warranty

88550 110v Knife
88550-BLADE Extra Heat Tip



Sports Mouthguard Trimming Bur

- Use to efficiently achieve smooth edges after trimming away excess plastic
- Follow with Polishing Brushes for fine contouring and smoothing
- 2.35mm shank; coarse grit

MGTB-1 Pack of 5



Finishing Liquid

- Rub on completed custom sports mouthguards for a high-quality finish
- Contains natural eucalyptus

4466N 250ml Bottle

Can only be shipped via ground delivery—sorry, no air delivery.

Achieve the perfect fit.

Appliance Trimming Kit

- Cordless Dremel 110v **DR001**
- Polishing Brushes Combo Pack **1898C**
- Essix Lab Knife **18800**
- 25 Refill Blades for Essix Lab Knife **188RF**
- 5 1/2" Mayo Scissors **18001**
- 3 1/4" Crown & Bridge Scissors **18002**

KIT-TRIM



Accentuator® Tool

- Helps Ensure a Snap-in Fit
- Use to guide plastic down into the interproximal undercuts

18943



Trim-Rite® Release Agent

- Helps avoid broken models
- Spray the cast to ensure a clean separation of the working cast from the thermoformed plastic
- Works with all thermoforming plastics
- Formulated not to turn plastic yellow

21004 12oz can

Can only be shipped via ground delivery—sorry, no air delivery.

Freeze Spray Coolant

- Provides enhanced adaptation
- Use to cool plastic after thermoforming to prevent malformation of plastic
- Without this rapid cooling, lingual or palatal edge may lift
- Also used to increase flexibility on thermal-activated arch wires
- Patent #6,371,759

22000 10oz can

22000-12 (12) 10oz cans

Can only be shipped via ground delivery—sorry, no air delivery.

Efficiently trim clear aligners and removable appliances.



Cordless Dremel

- Use to trim and polish removable appliances
- Variable speed - Seamlessly adjust to the ideal speed from 8,000 to 25,000 rpm
- Lithium Ion Battery - USB Charging with Battery Life Indicator
- Portable & Lightweight
- EZ Twist Nose cap for quick accessory changes
- Two-year manufacturer's warranty

D9001 Cordless Dremel 110v



Essix® Lab Knife

- Use for cutting or trimming plastic appliances

18800 Essix® Lab Knife

188RF 25 Refill Blades

18KIT Essix® Lab Knife & 25 Refill Blades



5 1/2" Mayo Scissors

- Ideal for cutting the labial outline on Essix® appliances with minimal jagged edges
- The curve of the scissors follows the curve of the arch
- Autoclavable

18001



4" Curve Scissors

- Help obtain a clean, smooth edge on the lingual anterior region of Essix® appliances
- Ideal for bleaching trays
- Autoclavable

18004



3 1/4" Crown & Bridge Scissors

- Designed for precise trimming, especially for the palatal or lingual on full arch appliances
- Autoclavable

18002



8" Universal Shears

- Excellent for quickly trimming excess plastic and more rigid materials
- Heavy-duty stainless steel
- Serrated edges
- Easy grip handles

18005

Help eliminate frayed edges.



Polishing Brushes

- Ideal for polishing edges of plastic appliances
- Scotch-Brite® material pre-mounted on mandrel
- Use with Cordless Dremel; max. speed 5,000rpm

- 18981** (12) Coarse
- 18982** (12) Medium
- 18983** (12) Fine
- 1898C** Combo Pack (4) of each

Scotch-Brite® is not a trademark of Dentsply Sirona.



Border Polishers

- Smooth edges of appliances
- Use with Cordless Dremel or slow-speed handpiece at a maximum of 10,000rpm

- 18973** Coarse Grit
- 18972** Medium Grit
- 18971** Fine Grit
- 18970** Set (1) ea. grit

Trimming Burs

- Trim plastic appliances easily and effectively
- Straight handpiece only; max. speed 10,000rpm

- 18979N** Narrow, Head Size: 2.3mm x 23mm
- 18979W** Wide, Head Size: 7mm x 14mm
- 18979S** Set of Wide & Narrow



Trimming Wheel

- Trim plastic appliances with this high-quality, rigid one-piece tool
- 10mm head size; 0.5mm thickness
- Max. speed 10,000rpm

18320



Micromotor

- Use with burs and polishing brushes to finish removable plastic appliances quickly
- 0-35,000rpm non-stage speed system
- Variable speed foot control; forward and reverse
- One-year warranty

207MM-110 110v/220v



Interproximal Reduction



THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

Power up your IPR.

SpaceFile® Contra-Angle Handpiece and Air Motor

- Use with ALL SpaceFile® Files
- Compact design makes reaching posterior regions easy
- Use with most standard slow-speed air motors and compressed air units such as A-dec, KaVo, W&H, NSK, and Midwest
- Reciprocating motion of 1.6mm
- Autoclavable
- One-Year Warranty

40027ANG SpaceFile Contra-Angle Handpiece

40027MOT SpaceFile Air Motor



SpaceFile® Interproximal Reduction Kit






- Precision Swiss-made, diamond-impregnated files accomplish a wide range of tasks including breaking contact, reduction, contouring, fine trimming, and finishing
- Can be used three ways, depending on your personal preference and case requirements
- Diamond-impregnated



Assorted Packs Includes 8 Files (1 of each thickness), 1 Serrated Contact Point Saw, and 1 Removable Grip	DOUBLE-SIDED	SINGLE-SIDED LEFT	SINGLE-SIDED RIGHT
	40034	40034-LS	40034-RS

Single Packs

Includes 8 files (all same thickness) and 1 Removable Grip

		THICKNESS		
Polishing A final step in tooth surface restoration.		40035-G	40035-LS-G	40035-RS-G
		0.085mm	0.070mm	0.070mm
Finishing For finishing prior to final polish.		40035-P	40035-LS-P	40035-RS-P
		0.110mm	0.100mm	0.100mm
XX Fine Pre-polish and finish reduced areas. Also use as initial strip for excessive crowding.		40035-W	40035-LS-W	40035-RS-W
		0.140mm	0.123mm	0.123mm
X Fine Pre-polish stripped areas. Also use as initial strip for crowded cases.		40035-Y	40035-LS-Y	40035-RS-Y
		0.160mm	0.138mm	0.138mm
Fine Use to initiate reduction in crowded cases. Also use to contour after reduction.		40035-R	40035-LS-R	40035-RS-R
		0.270mm	0.180mm	0.180mm
Medium Use for reduction with some crowding or when some space exists. Smaller grain can achieve results before using Coarse strip.		40035-BL	40035-LS-BL	40035-RS-BL
		0.330mm	0.210mm	0.210mm
Coarse Primary working strip for IPR. Also use to trim and contour.		40035-GR	40035-LS-GR	40035-RS-GR
		0.370mm	0.235mm	0.235mm
X Coarse Use for maximum reduction at start of IPR.		40035-BLK	40035-LS-BLK	40035-RS-BLK
		0.490mm	0.290mm	0.290mm
Serrated Contact Point Saw Use to break contact points on enamel, composite, amalgam, porcelain or gold.			40035-CS	0.300mm


SpaceFile® Removable Grip

40035-GRIP 3-Pack



Flexview® Interproximal Strips

- Optimal visibility
- Gateway area allows easy access for IPR stripping and finishing
- Autoclavable in aluminum tube **53500**



Item#	NARROW - 2.5MM WIDE / 148MM LONG			WIDE - 3.75MM WIDE / 148MM LONG				
	Single Packs 10 Strips (all same Grit)			Assorted Packs 15 Strips	Single Packs 10 Strips (all same Grit)			Assorted Packs 15 Strips
	5VNM1	5VNF1	5VNX1	5VN15	5VWM1	5VWF1	5VWX1	5VW15
Grit	Medium 45 micron	Fine 30 micron	Extra-Fine 15 micron	(5 of each Grit)	Medium 45 micron	Fine 30 micron	Extra-Fine 15 micron	(5 of each Grit)
Sided	Single	Single	Single		Single	Single	Single	
Thickness	0.13mm	0.10mm	0.08mm		0.13mm	0.10mm	0.08mm	

Diamond Interproximal Strips

- Gateway area allows easy access to contact points
- Stainless steel with electronically bonded diamonds
- Autoclavable in aluminum tube **53500**



Item#	NARROW - 2.5MM WIDE / 147MM LONG			WIDE - 3.75MM WIDE / 147MM LONG		
	Single Packs 10 Strips (all same Grit)			Single Packs 10 Strips (all same Grit)		
	53002	53001	53003	53101	53102	53103
Grit	Medium 45 micron	Medium 45 micron	Fine 30 micron	Medium 45 micron	Fine 30 micron	Fine 30 micron
Sided	Double	Single	Single	Double	Single	Single
Thickness	0.13mm	0.10mm	0.08mm	0.13mm	0.10mm	0.08mm



Qwikstrip™ Interproximal Strips

- Help Prevent Soft Tissue Irritation
- IPR strips allow IPR to be performed in the anterior and posterior of the mouth
- Serrated strip is used to break contact points prior to performing IPR
- Grip handles allow more control, so there is less chance of cutting patients' lips or gums
- Autoclavable

Grit	Single Packs 10 Strips (all same Grit)					Assorted Packs 10 Strips	Serrated Packs 10 Strips
	EXTRA COARSE 180 micron	COARSE 120-135 micron	MEDIUM 100-120 micron	FINE 45 micron	SUPERFINE 35 micron	(2 of each Single Pack Grit)	SERRATED
Single-Sided	QS-XCB	QS-CG	QS-MB	QS-FR	QS-SFY	QS-ASRT	QS-SW
Double-Sided	QSDS-XCB	QSDS-CG	QSDS-MB	QSDS-FR	QSDS-SFY	QSDS-ASRT	-
Curved Single-Sided	QSC-XCB	QSC-CG	QSC-MB	QSC-FR	QSC-SFY	QSC-ASRT	-

QwikStrip is not a trademark of Dentsply Sirona.

Autoclavable Aluminum Tube

- Sterilize up to 15 interproximal strips at a time
- 53500**



Serrated Interproximal Strips

- Use to open bonded proximal contact points
- Ultra-thin; made of stainless steel
- Thickness 0.05mm, width 6mm, length 150mm
- Autoclavable in aluminum tube 53500

53400 10 per pack



IPR Gauge Set

- Use to measure and confirm space
- Set of six stainless steel gauges
- 0.1mm, 0.2mm, 0.25mm, 0.3mm, 0.4mm, and 0.5mm
- Autoclavable
- Patent #6,413,086

IPR01





Posterior & Anterior ARS Stripping Burs

- Used for rapid enamel reduction
- Open field is required. Use separators first unless teeth are parallel
- Autoclavable; do not immerse in cold-sterilization solution




Friction Grip	Posterior ARS Cross-cut Carbide Stripping Burs	Anterior ARS Diamond Bur
Pack of 5	699C	55000
Tip Type	Safe	—
Head Size	1.0mm	1.6mm
Head Length	5.9mm	5mm
Max Speed	450,000rpm	<300,000rpm

Safe- And Standard-Tipped Burs

- Diamond coating helps produce a natural finish
- Autoclavable; do not immerse in cold-sterilization solution



Friction Grip	Safe-Tipped Finishing Burs			Standard-Tipped Finishing Burs			
Pack of 5	5MDST	5FDST	5EDST	848MD	514FD	504ED	848FD
Particle Size	Medium	Fine	Extra Fine	Medium	Fine	Extra Fine	Fine
Head Size	1.4mm	1.2mm		1.2mm	0.9mm		1.4mm
Head Length	10mm	8mm	10mm	10mm	3mm		10mm
Max Speed	<450,000rpm	300,000rpm		<450,000rpm	<40,000rpm	<30,000rpm	<40,000rpm



Standard-Tipped Bur Kit

Includes storage block and 12 burs:
 (2) Anterior Diamond ARS Burs (55000)
 (4) Posterior ARS Carbide Stripping Burs (699LC)
 (1) Medium Diamond Finishing Bur (848MD)
 (2) Fine Diamond Finishing Burs (514FD)
 (2) Extra-Fine Diamond Finishing Burs (504ED)
 (1) Extended Head Diamond Finishing Bur (848FD)
51200

Safe-Tipped Bur Kit

Includes storage block and 10 burs:
 (2) Anterior ARS Diamond Burs (55000)
 (4) Posterior ARS Carbide Stripping Burs (699C)
 (2) Fine Diamond Finishing Burs (5FDST)
 (2) Medium Diamond Finishing Burs (5MDST)
512C

- These easy-to-store kits have the burs needed for ARS - from start to finish. Individual burs may be autoclaved; do not immerse in cold-sterilization solution.



Patient Aids



Complete appliance cleaning solutions for your patients.



Retainer Brite® Cleaning Tablets

- Uniquely formulated to help preserve clarity and cleanliness
- It's easy to supply your patients
- Stock 1- or 3-month boxes to sell or give to patients
- Include a Patient Pak with a new appliance!

RB-36-CASE

(24) 1-Month Supply Boxes 36 tablets per Box

RB-92-CASE

(12) 3-Month Supply Boxes 96 tablets per Box

RB-TP

(100) Patient Paks. Each Pak has 2 tablets



Sonic Cleaners

- Makes dental appliances sparkling clean!
- Works great when used with Retainer-Brite
- Includes a six month replacement warranty
- 3.5" x 4", White
- Uses 2 AA batteries (included)

Non Personalized **RB-USC**
Personalized **RB-USC-C**



Dentsply Sirona provides high quality solutions for better, safer and faster dental care

Dentsply Sirona
Orthodontics

Chewies™ Aligner Tray Seaters

- Same size (13/8" l x 7/16" w) as standard dental cotton roll
- Made of soft plastic
- Patient Pak contains 2 seaters

AC-(color) Bag of 10 patient packs



Retainer Retrievers

- Remove appliances cleanly and efficiently
- Stores in appliance case
- 10 individual tamper-proof patient bags per pack


RR-1(color) Tropical Blue (TSB), Tropical Green (TSG)
 Tropical Orange (TSO), Tropical Pink (TSP)
 White (W)



Appliance Remover Tool

- Ideal for patients with large fingers
- Convenient; fits inside appliance cases, non-disruptive to existing protocol
- 40° Angle, 1 3/4 in, 3/4 in
- Patent- #7,011,517

ART01 Sold in Packs of 25



Therapy Wafers

- Helps reduce discomfort associated with wearing orthodontic braces
- Chewing on wafer may increase blood circulation around recently moved teeth by stimulating and massaging gum tissue

EB25 Cherry Scented, Clear
EB-25G Grape Scented, Purple
EB-25PBG Bubble Gum Scented, Pink
EB-25M Mint Scented, Green

Pack of 25





We're excited to share how you can put our high quality products and solutions to work for you to expand your services and grow your practice.

For more information contact your Dentsply Sirona Representative or visit us online at [dentsplysirona.com](https://www.dentsplysirona.com)

© 2021 Dentsply Sirona. All Rights Reserved. 0290-MKTMAT-522 Rev. 01 Issued 04-2021
Dentsply Sirona Orthodontics Inc • 7290 26th Court East • Sarasota, FL 34243

THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona